

V200 XL

Optimised settings

The V200 XL is the second model of the V-series in XL version. Made of lightweight aluminum to limit the total weight and therefore manageable for easy transport.



Heavy duty wheelchair

The **V200 XL** offers many technical customizations to personalize your **V200 XL** and perfectly fit your body and your wishes. With a maximum user weight of 170 kg, this V200 XL is one of our lighter wheelchairs for people with a higher body weight.

Colours

C59

- ▶ 3 Seat widths
- ▶ 4 seat heights
- ▶ Fixed seat depths
- ▶ 3 armrest positions
- ▶ Maximum user weight: 170 kg
- ▶ Different options

Technical data

Total width	740 / 770 / 810 mm
Total height	820 - 1000 mm
Total length	1050 mm
Folded width	385 mm
Seat width	530 / 560 / 600 mm
Height armrests	220 - 240 mm

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **V200 XL** complies with the highest norms, laid down by certified test institutes in most European countries.



Seat height	440 - 530 mm
Height backrest	400 mm
Seat depth	490 mm
Weight	20,15 kg
Maximum user weight	170 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you