


Easy Life





Jinny

All our shower benches are ideal for usage in the shower, bathtub, bathroom or other humid environments.

The specially designed rubbers at the end of each leg will create a safe stability for the bench and avoid the sliding of the bench.

width: 35 cm

depth: 35 cm

seat height: 33 cm - 51 cm

diameter seat: 32 cm

weight: 1,5 kg

users weight: 113 kg



Lilly



User-friendliness in combination with functional design.
All seats are equipped with a special anti-sliding top layer providing extra safety.
Available with or without a preshaped backrest.

width: 51 cm
depth: 51 cm
height: 77 - 87 cm
height back: 43 cm
seat width: 51 cm
seat height: 36 - 46 cm
weight: 3 kg
users weight: 113 kg

Jill



width: 51 cm
depth: 42 cm
seat height: 35 - 53 cm
seat width: 51 cm
weight: 2 kg
users weight: 113 kg



Melbourne

width: 48 cm
depth: 56 cm
height: 72 - 82 cm
height back: 37 cm
seat width: 41 cm
seat height: 34 - 44 cm
weight: 2,8 kg
users weight: 113 kg



Sidney

width: 50 cm
depth: 52 cm
seat height: 34 - 44 cm
seat width: 41 cm
weight: 2 kg
users weight: 113 kg



Kate & Katy

width: 80 cm
depth: 56 cm
height: 80 - 93 cm
height back: 42 cm
seat width: 66 cm
seat height: 41 - 54 cm
weight: 5,2 kg
users weight: 113 kg



width: 83 cm
depth: 61 cm
height: 80 - 91 cm
height back: 37 cm
seat width: 69 cm
seat height: 47 - 58 cm
weight: 4,5 kg
users weight: 113 kg





Ivy

Functional and sturdy bathtub rail providing extra support and safety whilst entering or exiting your bathtub.

Easy to install and applicable on most types of bathtubs due to its adjustable width.

width: 18 cm

depth: 23 cm

height: 48 cm

min. - max. opening: 6 - 11,5 cm

weight: 1,6 kg

users weight: 113 kg

Grace & Iris




The shower bench **Grace** has two handgrips and a back plate (to mount it to a wall). It can be folded away and locked into an upright position. Height adjustable legs are included when extra support is needed.

width: 48 cm
depth: 41 cm
height: 42 - 53 cm
seat width: 48 cm
seat height: 38 - 49 cm
weight: 3 kg
users weight: 113 kg

The foldable toilet bar **Iris** is designed to help you getting up from the toilet. It can be folded upwards if needed. Additionally, a toiletroll holder is provided.

width: 10 cm
length: 74 cm
height: 23 cm
weight: 1,8 kg
users weight: 113 kg



Liddy



The **Liddy** toilet rails provide additional comfort and support on your toilet. Adjustable in height and easy to install.

width: 56 cm
depth: 47 cm
height: 63 - 74 cm
weight: 1,45 kg
users weight: 113 kg

Tina & Stacy



The Vermeiren raised toilet seat **Tina** was designed to raise the height of a toilet to a comfortable level, making it easier to sit down or get up.

width: 52 cm
depth: 43 cm
height: 38 cm
seat width: 45 cm
width adjustment: 29 - 38 cm
weight: 2,5 kg
users weight: 136 kg



width: 57 cm
depth: 51 cm
seat width: 45 cm
seat height: 44 - 59 cm
height armrests: 23 cm
total height: 80 - 95 cm
weight: 8,6 kg
users weight: 113 kg





Izzie

The bathroom bars are available in different lengths from 30 to 80 cm.
With exception of the **Izzie R** and Izzie L, which are bars in L-shape.
They are designed as general support elements.
They are useable everywhere, including humid environments.

Izzie 30, 40, 45, 60, 80
length: 30 - 40 - 45 - 60 - 80 cm
depth: 7 cm
weight: 0,4 - 0,8 - 0,5
0,65 - 0,8 kg

Izzie L/R
length: 84 cm
height: 48 cm
depth: 7 cm
weight: 1,1 kg





Pepi



The **Pepi** bath deckchair for children was made in order to satisfy your expectations in the sector of products supporting children's hygiene.

The construction of the PEPI deckchair is made of aluminum pipes that ensure stable and failure-free use for a long time.

Thanks to its versatility of positioning, the product can be fully adapted to the user's needs.

width: 38 cm


seat depth: 20 cm

length of the total lying area: 110 cm

length of the support: 63 cm

max. users weight: 30 kg

weight: 3,7 kg



The helping hand **Nicky** is designed for daily assistance to grab or pick up small objects on the floor. The front of the helping hand Nicky is equipped with a small magnet. It can help you to pick up small magnetic objects such as keys, coins or cans.

Nicky

thickness: 1,2 cm

length: 76 cm

height: 11 cm

weight: 0,15 kg

colors: black, red, green, blue






Canes

Vermeiren offers a wide variety in canes for in- and outdoor usage. From stylish foldable models to fixed models.

Find the color and design that suits you most and have more enjoyment of your daily walks.



Medium: from 114 cm to 134 cm
Tall: from 134 cm to 154 cm

Ashley Tall & Medium (in height adjustable crutch)

- Beatrix 0,35 kg
- Betty 0,40 kg
- Bonny 0,35 kg
- Britney 0,35 kg
- Mona* 0,35 kg
- Mary* 0,40 kg
- Maud* 0,40 kg
- Monica* 0,35 kg



height: from 79 cm to 102 cm
(min. - max. height)
foldable or fixed model
users weight: 100 kg (* = 113 kg)



Beatrix Betty Bonny Britney Mary Maud Mona Monica
- foldable - - fixed -



Zoë

height: from 75 cm to 98 cm
base: cm. 26 x 26 x 30
weight: 0,60 kg

Rollator

A walker or rollator is ideal for people who require some additional support during their walks.

They are comfortable for in- and outdoor usage, to take a small rest during your walk or to carry something with you on the tray or in the shopping basket.

286B

width: 59 cm
height: 80 - 99 cm
depth: 70 cm
seat height: 61 cm
seat width: 43 cm
weight: 11 kg
users weight: 130 kg



275

width: 59 cm
height: 81 - 96 cm
depth: 76 cm
seat height: 58 cm
seat width: 44 cm
weight: 8,8 kg
users weight: 120 kg



275, one of the lightest models on the market

with elbow support

276



201 Delta

width: 63 cm
height: 87 - 99 cm
depth: 58 cm
weight: 7 kg
users weight: 130 kg



Streeter

Bright, colorful walkers with a functional design and ease of use. Completely foldable for transport in a car.

Available in different sizes according to your bodylength.

Mini

height: 46-57 cm
depth: 74 cm
folded: 20 cm
width basis: 53 cm
width top: 34 cm
weight: 4,75 kg
max. users weight: 50 kg

Medium

height: 47-60 cm
depth: 90 cm
folded: 17 cm
width basis: 55 cm
width top: 36 cm
weight: 5 kg
max. users weight: 65 kg

Large

height: 60-73 cm
depth: 96 cm
folded: 21 cm
width basis: 59 cm
width top: 40 cm
weight: 5,7 kg
max. users weight: 80 kg

X-large

height: 72-90 cm
depth: 102 cm
folded : 21 cm
width basis: 63 cm
width top: 43 cm
weight: 6,4 kg
max. users weight: 90 kg



Lyn & Lyna



Lyn

width: 58 cm
depth: 48 cm
folded: 11 cm
height: 81 - 99 cm
weight: 2,8 kg
users weight: 113 kg



Lyna (articulating)

width: 58 cm
depth: 48 cm
folded: 11 cm
height: 81 - 99 cm
weight: 2,8 kg
users weight: 113 kg

Options



Maureen



Amy



April

Lindsey
width: 63 cm
depth: 54 cm
folded: 12 cm
height: 79 - 97 cm
weight: 2,75 kg
users weight: 113 kg



Lindsey

Nightchair

Vermeiren offers a wide range in nightchairs, from fixed models to rolling models.

Soft and comfortable padding ensures a pleasant and natural sensation.



9062

width: 62 cm
depth: 61 cm
height: 85 cm
height armrests: 23 cm
height back: 40 cm
seat height: 50 cm
seat depth: 47 cm
seat width: 46 cm
weight: 8,7 kg
users weight: 120 kg

139B

width: 57 cm
depth: 82 cm
height: 95 cm
height armrests: 24 cm
height back: 41 cm
seat height: 52 cm
seat depth: 42 cm
seat width: 42 cm
weight: 12,9 kg
users weight: 100 kg



9060

width: 57 cm
depth: 56 cm
height: 82 cm
height armrests: 23 cm
height back: 36 cm
seat height: 50 cm
seat depth: 47 cm
seat width: 46 cm
weight: 7 kg
users weight: 120 kg





Eagle 620-625



The Vermeiren patientlift was designed to offer a real comfort for the patient itself and the nurse during the patients transfers.

Vermeiren designed this patientlift with specific attention for its functionality, safety and reliability.

Quality wheels of 100 mm ensure that the patientlift can be used on almost any surface.



Users weight


Eagle 620 : 150 kg

Eagle 625 : 175 kg

HYGIENIC MODEL

Available in
M - L - XL (solid & net)


	A	B	C	D
M	78	25	60	30
L	84	26	70	36
XL	91	27	82	44



COMFORT MODEL

Available in
M - L - XL (solid & net)


	A	B	C	D	E
M	78	40	45	30	39
L	84	45	50	36	41
XL	91	51	56	44	42



STANDARD MODEL

Available in
M - L - XL (solid & net)

	A	B	C	D
M	78	40	45	30
L	84	45	50	36
XL	91	51	56	44



Ramps

For a convenient “here-and-now” solution. Ramp kits that can be adapted to any situation.

All you need to do is to measure the height of the obstacle or the treshold and decide which type of ramp kit is suitable for you.



691



Model	Length	Width	Height	Weight	Max. weight
Kit nº.1	26 cm	75 cm	1 to 4 cm	1,3 kg	200 kg
Kit nº.2	51 cm	75 cm	4,5 to 7,5 cm	4,3 kg	200 kg
Kit nº.3	76 cm	75 cm	7,5 to 11,5 cm	9,2 kg	200 kg
Kit nº.4	101 cm	75 cm	11,5 to 15 cm	15,8 kg	200 kg
Telescopic ramps	106 - 190 cm	25 cm	5 cm	6,88 kg x 2	220 kg
	165 - 333 cm	25 cm	5 cm	10,6 kg x 2	220 kg

Kits 2,3 and 4 are delivered with an installation CD-rom.



V-drive

new

The V-Drive system simplifies your movements with a manual wheelchair. The V-Drive takes over all the work from your pushing assistant. Taking a steep hill or inclinations, bumps on the road or just a small journey outside... all this is just ideal for the V-Drive.

Vermeiren offers two versions of the V-Drive. A standard version (for people up to 135 kg) and a Heavy-Duty (H.D) version (for people up to 200 kg).

seat widths: 39 - 48 cm
(50 - 60 cm H.D.)
max. speed: 6 km/h
maximum inclination: 8°
action radius: ca. 15 cm
motor: 150 W (200 W H.D.)
batteries: 2 x 12 V 12 Ah
weight: 11,8 kg (w/o battery)
users weight: 135 kg (200 kg H.D.)



Geriatric chairs

NORMANDIE

Adjustable back/legrest
width: 66 cm
length: 75 cm
height: 120 cm
height back: 78 cm
seat width: 54 cm
seat height: 51 cm
weight: 26 kg
users weight: 130 kg



For optimal comfort at home, carecenter or hospital, these geriatric chairs have a wide variety of options.

Serumholders, tables, extra headsupport...

Make your choice from different colors and accessories.

NORMANDIE XXL

Adjustable back/legrest
width: 89 cm
length: 75 cm
height: 120 cm
height back: 78 cm
seat width: 74 cm
seat height: 51 cm
weight: 38 kg
users weight: 350 kg



LANGUEDOC

Adjustable back
width: 66 cm
length: 75 cm
height: 120 cm
height back: 78 cm
seat width: 54 cm
seat height: 51 cm
weight: 20 kg
users weight: 130 kg



PROVENCE

Fixed back
width: 66 cm
length: 75 cm
height: 120 cm
height back: 78 cm
seat width: 54 cm
seat height: 51 cm
weight: 19 kg
users weight: 130 kg



Alesia

Options:

Documentholder

Adjustable and removable headsupport

Legrest

Sidesupports

Table

Serum holder

Cannister holder

Urinebag holder



width: 66 cm
height: 133 cm
length: 78 cm
height armrests: 0 - 27 cm
height back: 80 cm
seat height: 53 cm
seat depth: 50 cm
seat width: 54 cm
inclination back: 30°
inclination seat: 12°
weight: 36,5 kg
users weight: 140 kg





ONTARIO

width: 81 cm
seat width: 45 cm
height armrests: 18 cm
seat height: 46 cm
total height: 112 cm
total length: 92 cm
seat depth: 50 cm
height back: 75 cm
weight: 57 kg
users weight: 200 kg

The Ontario liftchair offers you optimal comfort in your own house. With a simple push on a button you can relax, watch tv or have a little nap in your Ontario liftchair.



Artificial leather

Tissue



Liftchairs

The Vancouver and Montreal model will offer you the same possibilities as the Ontario. Make your choice out of different colors and textures for an optimal match with your existing interior.

MONTREAL

width: 79 cm
seat width: 45 cm
height armrests: 18 cm
seat height: 50 cm
total height: 109 cm
total length: 92 cm
seat depth: 53 cm
height back: 65 cm
weight: 58 kg
users weight: 150 kg

VANCOUVER

width: 73 cm
seat width: 45 cm
height armrests: 14 cm
seat height: 50 cm
total height: 109 cm
total length: 92 cm
seat depth: 50 cm
height back: 65 cm
weight: 56 kg
users weight: 150 kg

- ① The leg rest is operated while maintaining a straight back position.
- ② The back will fully recline in a relaxing position.



Artificial leather




Tissue



Beds

The Luna beds are stylish, high quality 4 section profiling electric beds with a user handset for convenience and full length wooden side rails for safety. The beds come complete with a lifting pole for extra user movement and independence. The multi adjustable 4 section mattress platform provides comfort and assists correct posture.



LUNA

width: 103 cm
length: 215 cm
min. height: 40 cm
max. height: 80 cm
inclination head: 83°
knee buckle: 16°
inclination feet: 30°
users weight: 140 kg

LUNA BASIC

width: 103 cm
length: 215 cm
min. height: 40 cm
max. height: 80 cm
inclination head: 83°
knee buckle: 16°
inclination feet: 30°
users weight: 140 kg



LUNA DELUXE

width: 103 cm
length: 215 cm
min. height: 40 cm
max. height: 80 cm
inclination head: 83°
knee buckle: 16°
inclination feet: 30°
users weight: 140 kg



LUNA ULTRALOW

width: 103 cm
length: 224 cm
min. height: 20 cm
max. height: 60 cm
inclination head: 83°
knee buckle: 16°
inclination feet: 30°
users weight: 140 kg

N.V. Vermeiren N.V.

Vermeirenplein 1-15

2920 Kalmthout

Belgium

Tel: 00 32 (0)3 620 20 20

Fax: 00 32 (0)3 666 48 94

info@vermeiren.be

Visit our website: www.vermeiren.com



Distributed by:

VERMEIREN

Opening up new horizons in mobility