



Detail A
Scale: 1:3

Drawing shows only right side

Notes:
*use only with 2993810 shortened

EP00266 / I

<h1>BACKREST TRIGO</h1>				Plan No.
				EP00266/I
Scale: 1:7	Drafted	Verified	Confirmed	
Date	06/11/2017	06/12/2023	06/12/2023	
Initials	AS	PM	AS	

Material (25-30)			
No.	Plan	Article	Quantity
3	YP05930	2785565	1
4	YP06265	2785265	1

Material (30-35)			
No.	Plan	Article	Quantity
3	YP05930	2785565	1
4	P00025189	2643065	1

Material (35-40)			
No.	Plan	Article	Quantity
3	YP05932	2786065	1
4	P00025189	2643065	1

Material (40-45)			
No.	Plan	Article	Quantity
3	YP05933	2786165	1
4	P00025189	2643065	1

Material (45-50)			
No.	Plan	Article	Quantity
3	YP05932	2786065	1
4	YP06269	2785465	1

Material				
No.	Plan	Article	Quantity	
5	YP05910 (25)	2786465	1	
	YP05911 (27)	2786565	1	
	YP05912 (29)	2786665	1	
	YP05913 (31)	2786765	1	
	YP05914 (33)	2786865	1	
	YP05915 (35)	2786965	1	
	YP05916 (37)	2787065	1	
	YP05917 (39)	2787165	1	
	YP05918 (41)	2787265	1	
	YP05919 (43)	2787365	1	
	YP05920 (45)	2787465	1	
	YP05921 (47)	2787565	1	
	YP05922 (49)	2787665	1	
	YP05906 (51)	2787765	1	
	6	YP05995 (25)	2789465	1
		YP05997 (27)	2789565	1
		YP05998 (29)	2789665	1
YP05999 (31)		2789765	1	
YP06000 (33)		2789865	1	
YP06001 (35)		2789965	1	
YP05992 (37)		2790165	1	
YP06002 (39)		2790065	1	
YP06003 (41)		2790265	1	
YP06004 (43)		2790365	1	

Material				
No.	Plan	Article	Quantity	
6	YP06005 (45)	2790465	1	
	YP06006 (47)	2790565	1	
	YP06007 (49)	2790665	1	
	YP06008 (51)	2790765	1	
	7	YP06896 (25)	2765765	1
		YP06897 (27)	2765865	1
		YP06898 (29)	2765965	1
		YP06899 (31)	2766065	1
		YP06900 (33)	2766165	1
YP06901 (35)	2766265	1		
YP06902 (37)	2766365	1		
YP06903 (39)	2766465	1		
YP06904 (41)	2766565	1		
8	YP06905 (43)	2766665	1	
	YP06906 (45)	2766765	1	
	YP06907 (47)	2766865	1	
	YP06908 (49)	2766965	1	
	YP06909 (51)	2767065	1	
	8	Fold. push handle	EP00246	1
	9	Y000462	F000050	1
	10	Y004353	F000411	1
	11	XP01560	2993561	4
12	X000273	G0606	1	

** When using a folded handle, select 2 sizes shorter tube.

EP00266 / I

BACKREST TRIGO			Plan No.	EP00266 / I
			Page 2/3	
Scale: 1:7	Drafted	Verified	Confirmed	
Date	06/11/2017	06/12/2023	06/12/2023	
Initials	AS	PM	AS	

Material			
No.	Plan	Article	Quantity
13	M6	2994414	2
14	M6x45	2994086	2
	M6x55*	2994376	2
15	---	---	---
16	---	---	---
17	CAP D16	1001598	1
18	M6x10	2994127	4
19	Hand grip GR-KS	2993809	1
	GR-ERGO-S	2993810	1
20	M6	2994063	1
21	YP06168(39)	2797065	1
	YP07800(41)	2782265	1
	YP07801(43)	2782365	1
	YP07802(45)	2782465	1
	YP07803(47)	2782565	1
	YP07804(49)	2782665	1
	YP07863(51)	2600665	1
22	Y017346	1748054	1
23	Y017412	1747800	1
24	Ø3x16 DIN1481	1003791	1
25	M6x15	1906622	1
26	X004876	17477	1
27	Handgrip	1906244	1

* Use with B82 fixed, B82S fixed and B82S C fixed

EP00266 / I

BACKREST TRIGO				Plan No.
Page 3/3				EP00266 / I
Scale: 1:7	Drafted	Verified	Confirmed	
Date	06/11/2017	06/12/2023	06/12/2023	
Initials	AS	PM	AS	