


Isometric view
Scale: 1:4

EP00351 / E


<h1>CROSS TRIGO 48-50</h1>				Plan No.
				EP00351 / E
Scale: 1:4	Drafted	Verified	Confirmed	
Date	18/10/2018	21/03/2022	21/03/2022	
Initials	PM	AS	TK	

Page 1/2

Material				
No.	Plan	Article	Quantity	
Width 48cm				
1	YP08546	2619265 (36cm)	1	
	YP08545	2619165 (38cm)	1	
	YP06707	2740165 (40cm)	1	
	YP06706	2740065 (42cm)	1	
	YP06705	2739965 (44cm)	1	
	YP06704	2739865 (46cm)	1	
	YP06703	2739765 (48cm)	1	
	YP05656	2720865 (50cm)	1	
	2	YP08548	2619465 (36cm)	1
		YP08547	2619365 (38cm)	1
YP06702		2739665 (40cm)	1	
YP06701		2739565 (42cm)	1	
YP06700		2739465 (44cm)	1	
YP06699		2739365 (46cm)	1	
YP06698		2739265 (48cm)	1	
YP05664		2721665 (50cm)	1	
Width 50cm				
1	YP08549	2619565 (36cm)	1	
	YP08550	2619665 (38cm)	1	
	YP08551	2619765 (40cm)	1	
	YP06697	2739165 (42cm)	1	
	YP06696	2738965 (44cm)	1	
	YP06695	2738765 (46cm)	1	
	YP06694	2738565 (48cm)	1	
	YP05657	2720965 (50cm)	1	
	2	YP08552	2620065 (36cm)	1
		YP08553	2619965 (38cm)	1
YP08554		2619865 (40cm)	1	
YP06693		2739065 (42cm)	1	
YP06692		2738865 (44cm)	1	
YP06691		2738665 (46cm)	1	
YP06690		2738465 (48cm)	1	
YP05665	2721765 (50cm)	1		

Material			
No.	Plan	Article	Quantity
3	YP05694	2720165	2
4	YP05996	27916	2
5	XP01851	2993570	1
6	XP01439	2993524	14
7	P00030394	2995366	4
8	X002136	1992414	4
9	M6	2994059	2
10	M8x115	2994008	1
11	NUTM8	2994009	1
12	M6 DIN125	2994068	2
13	ScrewD3.9X9.5	2994075	4
14	M6x45DIN912	2994058	1
15	M6x80DIN912	2994176	1

EP00351 / E

CROSS TRIGO 48-50				Plan No.
				EP00351 / E
Scale: 1:4	Drafted	Verified	Confirmed	
Date	18/10/2018	21/03/2022	21/03/2022	
Initials	PM	AS	TK	