



Isometric view  
Scale: 1:7

\*not depicted

Material 1995326 - Grey tyre			
No.	Plan	Article	Quantity
1	PUSH RIM	3993062	1
2	---	3993063	1
3	---	3993064	1
4	QR D15x131	3993065	1
5	QR D15x114	3993066	1
6	SPOKE	3990332	56
7	ANCHOR PLATE	3990348	8
8	SPACER	3990385	8
9	SCREW + WASHER	3992456	8
10	HANDGRIP	3992431	2
11	RIM 24"X1 3/8	3990339	2
12	PU TYRE 24 DG	3992910	2
13*	STRETCH FILM	3990378*	0,036
14*	WHEEL PLUG	1915211	2

Material 1917018 - Black tyre			
No.	Plan	Article	Quantity
1	PUSH RIM	3993062	1
2	---	3993590	1
3	---	3993591	1
4	QR D15x131	3993065	1
5	QR D15x114	3993066	1
6	SPOKE	3992911	56
7	ANCHOR PLATE	3992913	8
8	SPACER	3990385	8
9	SCREW + WASHER	3993159	8
10	HANDGRIP	3992431	2
11	RIM 24"X1 3/8	3992912	2
12	PU TYRE 24 BL	3992947	2
13*	STRETCH FILM	3990378	0,036
14*	WHEEL PLUG	1915211	2

EP00391 / B

HEM2-SET TYPE RSS B80

Plan No.

EP00391 / B

Scale: 1:7	Drafted	Verified	Confirmed	
Date	17/04/2019	29/10/2021	29/10/2021	
Initials	PM	DVB	TK	