



Material			
No.	Plan	Article	Quantity
2	XP01417 (A)	2706799	1
	X005632 (B)	2037299	1
3	XP01418 (A)	2706899	1
	X005633 (B)	2037399	1
4	Y022799	2048654	1
5	X001119	1001554	1
6	M6x16 DIN912	1001492	3
7	---	1005899	1
8	M8x30	1005860	1
9	M5x10 ISO7380	1007063	2
10*	Y020044	1790254	1
11*	D22	1990648	2
12*	M6x16 ISO7380	1008459	4

* For use on L34 / L14

(A) for tube $\Phi 22.2$

(B) for tube $\Phi 25.4$

* no. 10/11/12 do not use for Turios with L34/L14

E001090 / D

ATTENDANT CONTROLLER SE08 LINX				Plan No.
				E001090 / D
Scale: 2:5	Drafted	Verified	Confirmed	
Date	12/12/2017	29/06/2023	29/06/2023	
Initials	SD	AZ	PJ	