



RAIN COVER GEMINI II 32 2020: 2996835			
No.	Plan	Article	Quantity
1	NYLON BLACK	1008943	1,3m
2	PVC FOIL W:150CM	2996801	1,25m
3	LOGO VERMEIREN	3993014	1
4	POLYPROP BAND 215/25MM PP BLACK	1009516	0,5m
5	BORDERING TAPE 20MM BLACK	2993020	1
6	VELCRO BLACK 25 MM	1009709	0,01
7	VELCRO HOOK BLACK 25 MM	1009710	0,01

RAIN COVER GEMINI II 40 2020: 2996983			
No.	Plan	Article	Quantity
1	NYLON BLACK	1008943	1,5m
2	PVC FOIL W:150CM	2996801	1,45m
3	LOGO VERMEIREN	3993014	1
4	POLYPROP BAND 215/25MM PP BLACK	1009516	0,6m
5	BORDERING TAPE 20MM BLACK	2993020	1,4m
6	VELCRO BLACK 25 MM	1009709	0,02m
7	VELCRO HOOK BLACK 25 MM	1009710	0,02m

RAIN COVER GEMINI 2020				Plan No.
				EP00622 / A
Scale: 1:10	Drafted	Verified	Confirmed	
Date	29/09/2021	29/09/2021	29/09/2021	
Initials	OY	AZ	TK	

EP00622 / A