




Isometric view




Side view




Detail B



Top view



Detail C



Detail A

No.	Title	Art. nr.	Sticker	Quantity
1	Sticker Vermeiren 85X15	T001051		2
2	Sticker taxi hooks (SD 44-46: Pos.A)	T001076		4
	Sticker taxi hooks (SD 42: Pos.B)	T001076		4
3	Sticker serial number	T001097		1
4	Sticker name Eclips+	T001033		2

E000675 / D

STICKER CONFIGURATION ECLIPS+				Plan No. E000675 / D
STIKKER CONFIGURATIE ECLIPS+				
Scale: 3:40	Drafted	Verified	Confirmed	
Date	11/01/2011	27/01/2015	27/01/2015	
Initials	TT	TT	PJ	